




3 LAYER SURGICAL MASK - TIE (Type IIR)



- 3 Layer Surgical mask Type IIR
- Polypropylene (PP) Covered Soft / Adjustable Nose
- Adjustment with Tie
- CE mark, ISO, PLP



CODE	DESCRIPTION	PACKAGE
AMET2040	IIR SURGICAL FACEMASK	50 units/pack

DIRECTIVE

93/42/ECC Medical Device Regulation (Class II)

STANDARDS

EN ISO 13485:2016 Medical devices - Quality management systems - Requirements for regulatory purposes

EN ISO 9001:2000 Quality management systems. Requirements

EN ISO 10993-1 Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing

EN 14683:2005 Surgical masks - Requirements and test methods

EN 15223-1 Medical Devices – Symbols to be used with Medical Device Labels, Labelling and Information to be supplied – Part 1: General Requirements

OPTIONS

COLORS *(more options available)*

